

A different problem occurs with DNA that contains poly-A regions and dinucleotide repeats. The repeating area of DNA means the DNA polymerase is no longer able to bind to the DNA, resulting in enzymic slippage. The result of a sample containing any of these problems is often considered a failure and is unusable by the customer.

We experimented using samples of DNA from various customers, all of which we believed to either be GC rich or suffering from a poly-A region or a dinucleotide repeat. All of the samples with suspected secondary structures showed a significant improvement when repeated with SSS and produced strong reads with high quality well-resolved base peaks and with read lengths of approximately 800bp (figures 2c and d). The samples suspected of containing Poly-A or dinucleotide repeats showed no significant improvement when using the SSS reagent.

To summarise, we suggest that it is straightforward to resolve sequences where secondary structure has previously proved obstructive. We have now introduced the SSS reagent into the Geneservice Ltd portfolio and can cost effectively sequence these loci. For further details on sequencing service provision please go to:

www.geneservice.co.uk/service/sequencing/

Reference

Tarpinian S. Using the TripleMaster™ PCR system for Robust Amplification of GC Rich DNA Templates.

<http://www.eppendorf.com/script/cms-newspic.php?id=4479&inline=1&col=DOWNLOADFILE>




Figure 2. a) Raw control data using sample A and sequenced with ABI Big Dye v3.1, b) zoomed in electropherogram of control sequence, shows point at which secondary structure causes truncation of the sequence read c) Sample A sequenced using SSS reagent, note raw data no longer shows a truncated read d) zoomed in electropherogram of sample A sequenced with SSS reagent at same locus as in 2b, arrow indicates point at which truncation had previously occurred.